




Anti-Ticks and Fleas Device



GENERAL INFORMATION

Electronic device with ultrasonic emission which, using a small silicone support, can be applied to any collar or harness for dogs of all sizes. The ultrasonic emission, without any release of chemical substances, acts on the nervous system of ticks and fleas, keeping them away from the animal.

The ultrasounds emitted have a **range of 1.5 m** and are beyond the hearing ability of animals, the device is therefore totally harmless and does not cause any type of disturbance. The continuous operation of ultrasonic emission carries out a constant prophylactic action. Particularly suitable for protecting elderly or sick animals, puppies and pregnant females.



Its total harmlessness makes it particularly suitable for animals that live in close contact with adults and children.

TECHNICAL DATA

- Emission frequency of **40 KHz**;
- The device is tested at the **University of Camerino** which declares its effectiveness;
- **CE** tests for **ROHS** electromagnetic compatibility and ultrasonic emission at **TÜV** laboratories;
- **Waterproof** tests with **IPX4-IPX7** certificate;

- **CR2032** lithium battery with operating life of approximately **8/10 months in continuous use**;
- **Weight:** 20 g (0,7 oz)
- **Dimensions:** 40x45x14 mm (1,6x1,8x0,5 inches)

HOW TO USE

- Press  on the lid and release, 5 rapid flashes indicate that the product is working;
- During operation, the device automatically emits 2 flashes every 8 seconds to indicate correct functionality;
- To turn off the device, press  on the lid and release, a long flash indicates that the product is turned off;
- Battery replacement: Unscrew the three screws on the back of the device, once opened, remove the exhausted battery and replace with a new one.

IMPORTANT: during the parasite season, always leave the device running so as to obtain maximum effectiveness and better prophylaxis.



Anti-Ticks and Fleas Device

REPELLENT ACTIVITY OF ULTRASOUND AGAINST TICKS AND FLEAS

Long experiments at universities and in the field have demonstrated the effectiveness of ultrasound in repelling ticks and fleas from dogs and cats without the use of harmful and polluting chemicals. It should be underlined that the frequency used is around 40Khz! This means that it is absolutely not audible to dogs and cats who perceive frequencies up to 25Khz.

HOW DO THESE ULTRASOUNDS WORK?

In practice they disturb the nervous system of the parasites which lose their orientation and the strength to remain attached to the animal. The consequence is the gradual abandonment of the victim but above all the continuous prophylactic action as the weakened parasites move away from the animal and do not attack it.

In addition to its effectiveness, university tests demonstrate the total harmlessness of the product and the total absence of side effects. This means that the use of this product is important for dogs of all sizes and ages, sick animals, puppies and pregnant females.